

**POCT - TRAINING AND COMPETENCY RECORD**

Name of Trainee : \_\_\_\_\_  
 Designation : \_\_\_\_\_  
 Contact No. : \_\_\_\_\_  
 Location : \_\_\_\_\_

*N.B: This record should be kept by the Supervisor / Training Officer*

| No.             | Scope  | **Score | Date of Training | Signature of Trainer |
|-----------------|--|---------|------------------|----------------------|
| <b>1.0</b>      | <b>THEORY</b>  |         |                  |                      |
| 1.1             | Method/principle of test                                 |         |                  |                      |
| 1.2             | Handling and preparation of reagent/test strips/cassette |         |                  |                      |
| 1.3             | Proper storage conditions for strips/cassette/reagents   |         |                  |                      |
| 1.4             | Calibration  |         |                  |                      |
| 1.5             | Internal Quality Controls                                |         |                  |                      |
| 1.6             | External Quality Programme/Peer group comparison         |         |                  |                      |
| 1.7             | Maintenance  |         |                  |                      |
| 1.8             | Understanding error messages                             |         |                  |                      |
| 1.9             | Basic troubleshooting                                    |         |                  |                      |
| 1.10            | Test limitation(s) and interference(s)                   |         |                  |                      |
| 1.11            | Health and safety regulations & compliance               |         |                  |                      |
| 1.12            | Patient preparation requirements                         |         |                  |                      |
| 1.13            | Sample collection procedures                             |         |                  |                      |
| 1.14            | Result interpretation & clinical decision limits         |         |                  |                      |
| 1.15            | Incident reporting and adverse events                    |         |                  |                      |
| <b>2.0</b>      | <b>PRACTICAL SESSION</b>                                 |         |                  |                      |
| 2.1             | Maintenance  |         |                  |                      |
| 2.2             | Calibration  |         |                  |                      |
| 2.3             | IQC  |         |                  |                      |
| 2.4             | Sample collection procedures                             |         |                  |                      |
| 2.5             | Sample analysis  |         |                  |                      |
| TOTAL SCORE (%) |  |         |                  |                      |

| <b>**SCORE:</b>   | <b>TOTAL SCORE</b>  |
|-------------------|---|
| 1 = Poor          | >90 % = Competent, able to train and audit                        |
| 2 = Below average | 80 -90% = Competent and able to run test                          |
| 3 = Average       | 60 - 79% = Require supervision and re-training                    |
| 4 = Good          | <60% = Not competent and cannot perform test. Require re-training |
| 5 = Excellent     |   |

### SUMMARY / COMMENT

---



---



---

Assessed by :

Name :

Designation :

Date :

Verified by :

Name :

Designation :

Date :